

# F85: Optik Grundpraktikum

Jonathan Böcker, Vincent Mader

08.07.-10.07.2019

Heidelberg University

Tutor: Keerthan Subramanian

## **Abstract**

This FP lab course serves as an introduction to three important properties of light waves: Polarization, phase and frequency are observed using commonly used optical components such as wave-plates, polarizing beam-splitters, electro-optical modulators and acousto-optical modulators.

The material we use as an optical medium in this course is  $\text{LiNbO}_3$ .

# Contents

<b>1</b>	<b>Introduction</b>	<b>5</b>
<b>2</b>	<b>Theory</b>	<b>6</b>
2.1	Polarization . . . . .	6
2.1.1	Snell's law . . . . .	6
2.1.2	Malus law . . . . .	6
2.1.3	Brewster's angle . . . . .	6
2.1.4	Fresnel equations . . . . .	7
2.1.5	Reflection and polarization . . . . .	7
2.1.6	Birefringence . . . . .	8
2.1.7	Wave plates . . . . .	8
2.2	Electro-optical effect . . . . .	9
2.2.1	Pockels effect (linear electro-optical effect) . . . . .	9
2.2.2	Pockels effect in a non-isotropic crystal . . . . .	9
2.2.3	Pockels cell . . . . .	10
2.2.4	Faraday effect . . . . .	10
2.2.5	Optical isolator . . . . .	10
2.3	Acousto-optical effect . . . . .	11
<b>3</b>	<b>Experiment</b>	<b>13</b>
3.1	Polarizers and wave plates . . . . .	13
3.1.1	Brewster's angle . . . . .	13
3.1.2	Second glass plate . . . . .	13
3.1.3	Polarizers . . . . .	14
3.1.4	Beam splitters . . . . .	14
3.1.5	Half-waveplates . . . . .	14
3.1.6	Quarter-waveplates . . . . .	14
3.1.7	Setup from fig. 5 . . . . .	14
3.1.8	Mirror between quarter-waveplates . . . . .	15
3.1.9	Reflection back into quarter-waveplate . . . . .	15
3.1.10	Reflection back into half-waveplate . . . . .	15
3.1.11	Mirror . . . . .	15
3.1.12	Lamp light through optical isolator . . . . .	15
3.2	Electro-optical effect . . . . .	16
3.2.1	Measurement of amplification . . . . .	16
3.2.2	Calibration of Pockels cell . . . . .	16
3.2.3	Mach-Zehnder interferometer . . . . .	16
3.2.4	Manipulation of polarization and intensity . . . . .	18
3.2.5	Linear amplitude modulation . . . . .	18
3.3	Acousto-optical effect . . . . .	19
3.3.1	Experiments with a single AOM . . . . .	19
3.3.2	Experiments with two perpendicular AOMs . . . . .	20
<b>4</b>	<b>Results</b>	<b>22</b>
<b>5</b>	<b>Analysis</b>	<b>23</b>

<b>6</b>	<b>Summary</b>	<b>24</b>
<b>7</b>	<b>Discussion</b>	<b>25</b>

## List of Figures

1	refraction and reflection . . . . .	7
2	effect of a half-wave plate on polarization . . . . .	8
3	optical diode . . . . .	11
4	acousto-optical effect . . . . .	11
5	setups for waveplate experiment . . . . .	15
6	setup of Mach-Zehnder interferometer . . . . .	17
7	manipulation of polarization and intensity . . . . .	18

# 1 Introduction

## 2 Theory

### 2.1 Polarization

A beam of light can be described as a plane electromagnetic wave, like so:

$$\vec{E}(z, t) = E_0 \cdot \hat{e}_x \cdot e^{i(\omega t - kz)} \quad (1)$$

$$\vec{B}(z, t) = B_0 \cdot \hat{e}_y \cdot e^{i(\omega t - kz)} \quad (2)$$

The propagation vector  $\hat{e}_z$  together with  $\vec{E}$  and  $\vec{B}$  create an orthogonal system.

The direction of  $\vec{E}$  defines the polarization, which can be distinguished into three types:

1. linear polarization
2. elliptical polarization
3. circular polarization

In the case of elliptical and circular polarisation, the wave can be thought of as a superposition of two orthogonal plane waves with a phase difference  $\Delta\varphi = \frac{\pi}{2}$ : **What happens for  $\Delta\varphi \neq \frac{\pi}{2}$ ?**

$$\vec{E}_{\pm} = (E_x \hat{e}_x \mp i E_y \hat{e}_y) \cdot e^{i(\omega t - kx)} \quad (3)$$

If, in addition, both waves carry the same amplitude ( $E_x = E_y$ ), we speak of circular polarization.

#### 2.1.1 Snell's law

If a beam of light reaches the surface between two optical media with refractive indices  $n_1$  and  $n_2$ , the beam splits into a reflected and a refracted one. The refraction angle  $\alpha_2$  can be determined by

$$n_1 \cdot \sin(\alpha_1) = n_2 \cdot \sin(\alpha_2) \quad (4)$$

The angles  $\alpha_1$  and  $\alpha_2$  are taken between the incoming/outgoing beam of light and the surface.

For reflection, the condition  $\alpha_{in} = \alpha_{out}$  must hold.

#### 2.1.2 Malus law

The Malus law states that when a perfect polarizer is placed in a polarized beam of light, the intensities before and after transversing the polarizer are related by

$$I_f = I_i \cdot \cos^2(\theta) \quad (5)$$

Here,  $\theta$  is the angle between the light's initial polarization direction and the axis of the polarizer.

#### 2.1.3 Brewster's angle

Brewster's angle is a specific angle of incidence at which light with a particular polarization is perfectly transmitted through a transparent dielectric surface, without any reflection.

When unpolarized light is incident at this angle, the light that is reflected from the surface is therefore perfectly polarized. **Meaning? Sources?**

### 2.1.4 Fresnel equations

The Fresnel equations describe the reflection and transmission of light when incident on an interface between different optical media. **equations**

### 2.1.5 Reflection and polarization

Consider again a beam of light traversing from one medium into another. The beam splits into a reflected and a refracted part. From Maxwell's equation the amplitudes and intensities of the two beams can be derived. The amplitude coefficients corresponding to the reflected and the transmitted beams are labeled  $r$  and  $t$ . The intensity coefficients  $R$  and  $T$  can easily be calculated by taking the square of  $r$  and  $t$ . If the polarization is perpendicular to the plane of incidence, the effect is called transversal-electric polarization (German: S-Polarisation).

The coefficients are given by

$$r_{TE} = -\frac{\sin(\alpha_1 - \alpha_2)}{\sin(\alpha_1 + \alpha_2)} \quad (6)$$

and

$$t_{TE} = \frac{2 \cdot \sin(\alpha_1) \cdot \cos(\alpha_2)}{\sin(\alpha_1 + \alpha_2)} \quad (7)$$

If the wave is polarized parallel to the plane of incidence, we speak of transversal-magnetic polarization (German: P-Polarisation), the coefficients are given by

$$r_{TM} = \frac{\tan(\alpha_1 - \alpha_2)}{\tan(\alpha_1 + \alpha_2)} \quad (8)$$

and

$$t_{TM} = \frac{2 \cdot \sin(\alpha_1) \cdot \cos(\alpha_2)}{\sin(\alpha_1 + \alpha_2) \cdot \cos(\alpha_1 - \alpha_2)} \quad (9)$$




FIGURE 1.8

Figure 1: refraction and reflection

**plot?**

### 2.1.6 Birefringence

If the medium is perfectly homogenous and isotropic, the refraction index is the same in all directions. Here, there is the simple relation  $\vec{D} = \varepsilon \vec{E}$  with permittivity  $\varepsilon = n^2$ . In general, a medium is likely to be anisotropic. The permittivity is best described by a tensor  $\varepsilon_{ij}$ .

uniaxial, biaxial media

### 2.1.7 Wave plates

Wave plates are birefringent crystals used to shift the phase of a traversing wave by a certain fraction of a whole period.

Quarter-wave plates can be used to turn a light beam's polarization from elliptical to linear and the other way round. For this its width ideally is  $(m + \frac{1}{4}) \cdot \lambda, m \in \mathbb{N}$ .

Half-wave plates can be used to change the polarization state of linearly polarized light. The optimal width for this is  $(m + \frac{1}{2}) \cdot \lambda$ .




Figure 2: effect of a half-wave plate on polarization

## 2.2 Electro-optical effect

Some materials experience a change in their optical properties when they are brought into an external electric field. If the refractive index  $n$  is a function of the applied field  $E$ , we speak of an electro-optical modulator (EOM).

### 2.2.1 Pockels effect (linear electro-optical effect)

The refractive index  $n = n(E)$  can be expanded around  $E = 0$  for small field strengths.

With  $r = -\frac{2}{n^3} \left( \frac{dn}{dE} \right) \Big|_{E=0}$  and  $s = -\frac{1}{n^3} \left( \frac{d^2n}{dE^2} \right) \Big|_{E=0}$ , this leads to the relation

$$n(E) = n_0 - \frac{1}{2} r n^3 E - \frac{1}{2} s n^3 E^2 \quad (10)$$

The linear electro-optical effect, also called Pockels effect, occurs for  $r \gg s$ :

$$n(E) \approx n_0 - \frac{1}{2} r n^3 E \quad (11)$$

Here,  $r$  is called the Pockels coefficient. If, on the other hand,  $r \ll s$ , the quadratic dependence of  $n$  on  $E$  is known as the Kerr effect. This effect will not be studied in this lab course.

### 2.2.2 Pockels effect in a non-isotropic crystal

Since the electro-optical crystals are in general birefringent, the Pockels coefficient is not a scalar, but a tensor. For the material  $\text{LiNbO}_3$ , the entries of this coefficient-tensor are given by

$$\begin{pmatrix} r_{11} & r_{12} & r_{13} \\ r_{21} & r_{22} & r_{23} \\ r_{31} & r_{32} & r_{33} \\ r_{41} & r_{42} & r_{43} \\ r_{51} & r_{52} & r_{53} \\ r_{61} & r_{62} & r_{63} \end{pmatrix} = \begin{pmatrix} 0 & -r_{22} & r_{13} \\ 0 & r_{22} & r_{13} \\ 0 & 0 & r_{33} \\ 0 & r_{51} & 0 \\ r_{51} & 0 & 0 \\ -r_{22} & 0 & 0 \end{pmatrix} = \begin{pmatrix} 0 & -3.4 & 8.6 \\ 0 & 3.4 & 8.6 \\ 0 & 0 & 30.8 \\ 0 & 28 & 0 \\ 28 & 0 & 0 \\ -3.4 & 0 & 0 \end{pmatrix} \cdot 10^{-12} \frac{\text{m}}{\text{V}}$$

If no external field is applied to the crystal, the indicatrix is given by

$$\frac{1}{n_o^2} x^2 + \frac{1}{n_o^2} y^2 + \frac{1}{n_e^2} z^2 = 1 \quad (12)$$

In this lab course, the electric field is applied along the extraordinary axis of the Pockels cell, the light propagates along one of the ordinary axes. With  $\vec{E} = E \cdot \hat{e}_z$  the indicatrix can be written as

$$\left( \frac{1}{n_o^2} + r_{13} E_z \right) x^2 + \left( \frac{1}{n_o^2} + r_{13} E_z \right) y^2 + \left( \frac{1}{n_e^2} + r_{33} E_z \right) z^2 \quad (13)$$

For small applied field strengths this leads to the relations

$$n'_o(E) \approx n_o - \frac{1}{2} r_{13} n_o^3 E_z \quad (14)$$

$$n'_e(E) \approx n_e - \frac{1}{2} r_{33} n_e^3 E_z \quad (15)$$



### 2.2.3 Pockels cell

An electro-optical crystal between two capacitor plates is called an electro-optical modulator, abbreviated as EOM. A Pockels cell is an EOM which exhibits the Pockels effect, meaning that the refractive index is approximately proportional to the strength of an applied electric field  $\vec{E}$  (for small values of  $|\vec{E}|$ )

If light with a wavelength  $\lambda$  travels through a medium with refractive index  $n$  for a distance  $L$ , the accumulated phase is given by the relation

$$\Phi = 2\pi n \frac{L}{\lambda} \quad (16)$$

The phase difference between a beam of light traveling along the ordinary vs the extraordinary axis is thus:

$$\Delta\Phi = 2\pi(n_e - n_o) \frac{L}{\lambda} \quad (17)$$

With a Pockels cell, we can manipulate the magnitude of both  $n_o$  and  $n_e$ . Plugging equation 14 and 15 into this last equation, we get the following:

$$\Delta\Phi(E) = 2\pi \frac{L}{\lambda} \left( n_e - n_o - \frac{1}{2}(r_{33}n_e^3 - r_{13}n_e^3 - r_{13}n_o^3)E_z \right) \quad (18)$$

This can be written as

$$\Delta\Phi(V) = \Phi_0 - \pi \frac{V}{V_\pi} \quad (19)$$

with  $V = Ed$ ,  $\Phi_0 = 2\pi \frac{L}{\lambda}(n_e - n_o)$  and  $V_\pi = \frac{d}{L\lambda} r_{33}n_e^3 - r_{13}n_o^3$ . A Pockels cell can be used to modulate the intensity of light. For this, the cell has to be positioned between two crossed polarizers, each at an angle of  $45^\circ$  relative to the crystal's optical axis. The transmittance of this setup is

$$T(V) = \sin^2 \left( \frac{\Phi_0}{2} - \frac{\pi V}{2 V_\pi} \right) \quad (20)$$

### 2.2.4 Faraday effect

Some media become optically active when an axial magnetic field is applied. Optically active means that the plane of polarization of linear polarized light is rotated. The angle of rotation  $\alpha_{rot}$  depends on the length of the medium  $L$ , the magnetic field strength  $|\vec{B}|$  and the so-called Verdet constant  $v$ :

$$\alpha_{rot} = vL|\vec{B}| \quad (21)$$

The Verdet constant is a function of the wavelength  $\lambda$ :

$$v = -\frac{\pi\gamma}{\lambda n} \quad (22)$$

Here,  $\gamma$  is a material constant of the medium, the so-called magnetogyration coefficient.

### 2.2.5 Optical isolator

The Faraday effect can be used to build an optical isolator, also called an optical diode.




Figure 3: optical diode

## 2.3 Acousto-optical effect

If a sound wave passes through a crystal, its density varies periodically, which also leads to a periodic variation in the refractive index. A plane sound wave with wavelength  $\lambda_s$  in a crystal with initial refractive index  $n_0$  can be described by

$$n(x, t) = n_0 - \Delta n \cdot \cos\left(\omega t - \frac{2\pi}{\lambda_s} x\right)$$

Here, the amplitude  $\Delta n = \frac{1}{2}pn^3s_0$  depends on the photo-elastic constant  $p$  and the amplitude of the strain  $s_0$ . In this course the main interaction between sound and light will be the so-called Debye-Sears effect which occurs for short interaction lengths, i.e. for a thin crystal or thin sound beam. Parts of the light beam that are travelling through the denser regions experience a phase shift. Maxima in the far field can be observed, if these parts of the beam interfere constructively with each other, as can be seen in the next figure.




Figure 4: acousto-optical effect

Because the light is interacting with a moving sound wave, the diffracted light is Doppler-shifted:

$$\omega_{out} = \omega_{in} + m\Delta\omega = \omega_{in} + m\omega_s$$

A different approach to explain the effect of an AOM on a light beam is to describe the system as two scattering quasi-particles, namely a phonon and a photon. Their momentum vectors are  $\hbar\vec{k}_l$  and  $\hbar\vec{k}_s$ , respectively. Conservation of momentum and energy leads to the relations

$$\vec{k}_{l,f} = \vec{k}_{l,i} \pm m\vec{k}_s$$

and

$$\nu_{l,f} = \nu_{l,i} \pm m\nu_s$$

Here,  $m$  is the diffraction order, i.e. the number of phonons that interacted with the photon. Constructive interference occurs when

$$\sin\theta_1 + \sin\theta_2 = m\frac{\lambda}{\lambda_s}$$

Here,  $\theta_1$  is the angle under which light enters the crystal and  $\theta_2$  the diffraction angle. **efficiency**

## 3 Experiment

### 3.1 Polarizers and wave plates

Date of carrying out experiment: 8th of July, 2019

This first part of the experiment has the purpose of getting used to the following optical devices: polarizers, beam splitters and waveplates. The laser used in this part of the experiment has a wavelength of 632.8 nm and a power output of less than 1 mW.

#### 3.1.1 Brewster's angle

First, we would like to find Brewster's angle. This allows us to later on calibrate the optical devices. For this, we take a glass plate and mount it in front of the laser at an angle  $\alpha$ . This angle can now be varied. If the angle between the incoming beam and the surface of the glass plate equals Brewster's angle, then the reflected beam contains only photons that are polarized perpendicular to the plane of incidence.

This can in principle be verified easily by mounting a polarizer in the path of the reflected beam, but the polarizers are not yet calibrated. This means that we have to turn the polarizers by  $360^\circ$  for each angle  $\alpha$ . We're looking for the setup which lets us filter out the maximum amount of light. A minimum in intensity tells us that the polarizer filters out all or almost all of the light, which means the beam contains only a single polarization component. This is exactly what one would expect at Brewster's angle.

Unfortunately, it's not easy to do this in a precise way, since two quantities have to be varied at the same time: the angle  $\alpha$  and the angle of the polarizer. Also, the laser beam is polarized linearly. This means, that both the reflected and the refracted/transmitted beam are polarized. It is always possible to block the reflected ray using a polarizer. Because of this, we mount a quarter-waveplate between laser and glass plate to make sure both components are present.

Our result:

$$\alpha_{Brewster} \approx (60 \pm 10)^\circ$$

Now we want to observe what happens when we use other optical devices after the beam has been reflected at Brewster's angle. As expected, the beam splitter splits the beam into one beam with very high intensity and another beam whose intensity is almost zero.

A polarizer either blocks the beam or has no effect, depending on orientation.

#### 3.1.2 Second glass plate

A second glass plate is mounted perpendicular to the plane of incidence of the first glass plate, so that the beam reflected on the first glass plate is reflected off of it. The incoming beam is once again split into two. In contrast to the previous observation, it is now possible to use a polarizer to block both beams for any orientation of the glass plate. This is to be expected, since the beam that was reflected on plate 1 is polarized, so both beams after plate 2 are polarized as well.

### 3.1.3 Polarizers

Now we want to calibrate the polarizers, i.e. we want to know for which angles linearly polarized light is blocked. For this, we use a similar setup as in the last observation. The laser is reflected off the glass plate at Brewster's angle through the polarizer. By varying the orientation of the polarizer and looking for a minimum of the intensity, we find:

blocking angle for polarizer 130000:  $\alpha_1 = (100 \pm 5)^\circ$

blocking angle for polarizer 130200:  $\alpha_2 = (101 \pm 5)^\circ$

Because we know that the beam reflected from the glass plate is polarized vertically, we can conclude that the polarizers are calibrated in such a way that the vertically polarized laser light can pass through the polarizer if it is set to  $\approx 10^\circ$ ,

### 3.1.4 Beam splitters

We want to find out the polarization of the two beams exiting a beam splitter. After mounting a polarizer behind the beam splitter and looking for which angles the beams can pass through, we observe that both beams show maximum intensity for an angle of  $0^\circ$  and minimum for  $90^\circ$ .

### 3.1.5 Half-waveplates

Now we'd like to observe how a half-waveplate or a quarter-waveplate influences the polarization state of the light beam.

First, we take a look at the half-waveplate. It is mounted between the glass plate and a polarizer. Minima are observed for the angles  $2^\circ$ ,  $91^\circ$ ,  $182^\circ$  and  $270^\circ$  with uncertainties of about  $5^\circ$ . These are the positions where the incident polarization is rotated by  $90^\circ$ .

### 3.1.6 Quarter-waveplates

Herfore, a quarter-waveplate is mounted between two polarizers. For waveplate a, the minima are observed at the angles  $8^\circ$ ,  $98^\circ$ ,  $188^\circ$  and  $279^\circ$ . For waveplate b, the minima are observed for the angles,  $59^\circ$ ,  $151^\circ$ ,  $241^\circ$  and  $330^\circ$  with an uncertainty of about  $5^\circ$

### 3.1.7 Setup from fig. 5

Next, we realize the experimental setup shown in figure 5a. Waveplate a is set to the angle  $8^\circ + 45^\circ = 53^\circ$  and waveplate b to  $59^\circ + 45^\circ = 104^\circ$ . We find minima for the polarizer angles  $180^\circ$  and  $360^\circ$ . If the second waveplate is turned by  $180^\circ$  around the axis on which it is mounted, we get minima at  $88^\circ$  and  $266^\circ$ .

To realize the setup from figure 5b, we need to rotate the initial polarization by  $45^\circ$ . This can be done with a half-waveplate that is set to an angle that corresponds to halfway from a minimum to the next maximum. That is about  $22.5^\circ$ . The waveplate angle for a is now  $59^\circ$ , the angle for b is  $8^\circ$ . The minima are at angles of  $25^\circ$  and  $205^\circ$ . After turning waveplate b on its axis, the minima are at roughly the same positions as before. This is two be expected, **because...**




Figure 5: setups for waveplate experiment

### 3.1.8 Mirror between quarter-waveplates

Now, a mirror is positioned between the two quarter-waveplates in a way so that the angle of reflection is bigger than  $45^\circ$  (total angle bigger than  $90^\circ$ ). After this is done, we do not observe any dependence of the intensity on the polarizer angle. We conclude that the light is circularly polarized. *maybe both quarter-waveplates were switched, wrong angles*

### 3.1.9 Reflection back into quarter-waveplate

The laser beam is lead through a quarter-waveplate and onto a mirror. The mirror is positioned perpendicular to the optical axis, so that the beam goes back through the quarter-waveplate. If the quarter-waveplate is set to an angle of  $8^\circ$ , the minima are found for the polarizer angles  $165^\circ$  and  $350^\circ$ . If it is instead set to  $53^\circ$ , the minima are at  $85^\circ$  and  $265^\circ$

### 3.1.10 Reflection back into half-waveplate

The same as above is done again, this time with a half-waveplate. At a waveplate angle of  $2^\circ$ , no minima could be found. Then we set the waveplate angle to  $47^\circ$ . Here, we did find minima at  $75^\circ$  and  $255^\circ$ .

### 3.1.11 Mirror

A mirror is set up in a way such that the laser beam is reflected from the glass onto the mirror in an angle bigger than  $45^\circ$ . Minima are detected for the polarizer angles of  $90^\circ$  and  $270^\circ$ . Now we do the same for perpendicular polarization. Now the minima are at  $186^\circ$  and  $6^\circ$ .

### 3.1.12 Lamp light through optical isolator

Looking through the optical isolator at the ceiling lamp, the lamp appears either orange or green, depending on the orientation of the isolator. *Why?*

## 3.2 Electro-optical effect

Date of carrying out experiment: 9th of July, 2019

The laser used in this part of the experiment produces light with a wavelength of  $\lambda = 635$  nm and outputs a power of less than 1 mW. The material in the Pockels cell is  $\text{LiNbO}_3$ , which is a uniaxial crystal with  $3m$ -symmetry ( $n_1 = n_2 = n_o, n_3 = n_e$ )

Pockels coefficients are very hard to measure, can be affected by crystal impurities

### 3.2.1 Measurement of amplification

First, we want to determine the amplification factor of the high voltage power supply. The DC power supply is connected to the amplifier. For several DC voltages  $U_0$  we record the corresponding high voltage value  $U_{amp}$ .

$U_0$ [V]	6.50	6.00	5.50	5.00	4.50	4.00	3.50	3.00	2.50	2.00	1.50	1.00	0.50
$U_{amp}$ [kV]	1.95	1.80	1.65	1.50	1.35	1.20	1.05	0.90	0.75	0.60	0.45	0.30	0.15

We take care not to let the amplified voltage get higher than 2.00 kV, this is the maximum voltage the Pockels cell can endure. The uncertainty on the display of the power supply devices is 0.005 V for  $U_0$  and 0.005 kV for  $U_{amp}$ . other error?

### 3.2.2 Calibration of Pockels cell

To determine the directions of the optical axes in the Pockels cell, it is placed between two polarizers. The polarizers are oriented perpendicular to each other, so that without the cell no light can pass through.

polarizer serial number	polarizer angle
130000	100°
130200	10°

The laser light passes first through the polarizer 130200, then through the Pockels cell and 130000.

The Pockels cell angle scale goes from  $-90^\circ$  to  $+90^\circ$  and the minima are observed for the Pockel cell angles  $10^\circ$  and  $-80^\circ$ . Minima are observed when the laser beam is parallel to either the fast or the slow axis of the crystal.

### 3.2.3 Mach-Zehnder interferometer

We would like to quantify the effect that the Pockels cell has on the phase of a laser beam. For this, a Mach-Zehnder interferometer is set up. The setup can be seen in figure 6.




Figure 6: setup of Mach-Zehnder interferometer

After the separated laser beams are joined back together, we want to make sure they are as close to being parallel as possible. A piece of paper is held into the beam directly behind the beam splitter. The mirrors are oriented in such a way that both beams hit the paper in a single spot. Then the distance between the piece of paper and the beam splitter is increased and the mirror positions once again adjusted.

A triangular waveform is chosen on the function generator. The frequency of this signal is  $\approx 3.5$  Hz !!! The generator is connected to the amplifier, the amplifier to the Pockels cell. The output of the function generator as well as the photo diode are connected to the oscilloscope. The oscilloscope output can be saved to .csv and used for later analysis.

This is done for the Pockels cell angles  $-80^\circ$ ,  $10^\circ$ ,  $-90^\circ$  and  $90^\circ$



### 3.2.4 Manipulation of polarization and intensity

The setup is shown in figure 7. The angle between the polarization axes of the filters and the Pockels cell need to be  $45^\circ$ , because...

The angle of the first filter (set to vertical polarization) is  $10^\circ$ , the second one  $100^\circ$ . For the Pockels cell, we know from the calibration that an angle of  $+45^\circ/-45^\circ$  relative to the vertical polarization axis corresponds to an actual angle setting of  $55^\circ/-35^\circ$  at the Pockels cell.

The transmitted power on the photodiode is recorded as a function of the applied high voltage by choosing again a triangular waveform on the the function generator. This is done for the angles  $-45^\circ$  and  $45^\circ$ .




Figure 7: manipulation of polarization and intensity

### 3.2.5 Linear amplitude modulation

Now, the Pockels cell is connected to a DC high voltage. The amplifier is set to intern. Multiple measurements will be made, the DC voltage will be varied from one measurement to the next. A modulated signal is added to the DC voltage and the resulting signal is given to the Pockels cell. The frequency of the modulated (sinoid) signal is about  $1.0000 \pm 0.0005$  kHz.

For the following DC voltages and Pockels cell angles the waveform is saved:

Pockels cell angle	$-35^\circ$	$55^\circ$	$-90^\circ$	$90^\circ$
DC voltage	0.60 kV	0.46 kV	0.49 kV	0.45 kV

### 3.3 Acousto-optical effect

Date of carrying out experiment: 10th of July, 2019

#### 3.3.1 Experiments with a single AOM

A high frequency acoustic signal is applied to the crystal by connecting a DC voltage to the VCO.

A fairly symmetric interference pattern showing 5 maxima can be realized with the frequency-controlling voltage being set to 9.4 V and the VCO level-controlling voltage to 1.9 V.

For different values of the VCO frequency-controlling voltage  $U_f$ , we record the distance  $a_i$  between the maxima of order 0 and  $i$ . The AOM frequency can be calculated from the applied voltage. Since the range of the AOM is  $(110 \pm 25)$  MHz<sup>1</sup> and the voltage can be varied by  $\pm 10$  V, we conclude that varying the voltage by 1 V corresponds to a frequency change of 2.5 Hz.

The distance between the crystal and the screen is  $d = (194 \pm 3)$ .

$U_f$ [V]	1	2	3	4	5	6	7	8	9	10
$a_1$ [cm]	4.4	4.3	4.1	3.9	3.7	3.5	3.2	3.1	2.9	2.7
$a_2$ [cm]	-	-	-	7.8	7.4	7.1	6.6	6.4	5.9	5.6

The uncertainties of the positions are estimated to be about 2 mm, since the points in the interference pattern have a certain spatial extension.

Now, we want to compare the intensity of the 0th order maximum with the intensity of the 1st order maximum. For that, we mount a photodiode in the path of the laser beam. The photodiode voltage for the undiffracted beam is read from the oscilloscope as  $U_{0,1} = (96 \pm 2)$  mV if the diode is placed directly behind the laser or  $U_{0,2} = 89 \pm 2$  mV if the diode is placed at a certain distance from the cell ( $\approx 90$  cm). The distance between the cell and the laser is  $\approx 20$  cm.

$U_f$ [V]	1	2	3	4	5	6	7	8	9	10
$U_{max,1}$	46	49	$43 \pm 2$	49	77	$77 \pm 2$	77	70	$55 \pm 2$	63

The uncertainties for the voltages are  $\pm 1$  mV each, if not otherwise specified.

für große  $U_f$ : Fehler: Maxima werden Striche statt Punkte, nicht alles auf Diode

Because another group was in the room during the measurement, we couldn't turn off the lights. Because of this, a background measurement is done at the end. The photodiode voltage with the lights off is  $100 \pm 1$  mV.

Now, we want to do a similar measurement, but this time vary the amplitude  $A$  instead of the frequency  $f$  of the sound wave. This we do by varying the voltage  $U_A$  for the VCO level. The voltage controlling the frequency is set to  $U_f = (7.00 \pm 0.05)$  V.

$U_A$ [V]	0.20	0.40	0.60	0.80	1.00	1.20	1.40	1.60	1.80	2.00
$U_{max,1}$	-	$19 \pm 2$	$34 \pm 2$	$69 \pm 2$	79	83	82	81	83	85

The uncertainty of the voltage  $U_A$  is estimated to be about 0.05 V and stems from the reading process on the display, which can only show one digit after the comma. As before, the uncertainty of  $U_A$  is  $\pm 1$  if not otherwise specified.

nicht symmetrisch ausgerichtet: Problem?

### 3.3.2 Experiments with two perpendicular AOMs

For the next experiment, both AOMs are placed near each other ( $\approx 1$  cm) and perpendicular to each other. Each of them is connected over the amplifier to one of the outputs of the function generator with two outputs. The setup can be seen in figure ???. We expect to now see a two-dimensional interference pattern, which is also what we see.

If the offset voltages are increased, the fringes in the pattern move closer to each other in the corresponding direction (left/right or up/down depending in which AOM is controlled).

When we modulate the sound amplitudes, we see that the interference pattern starts blinking. This is to be expected, since the amplitude of the sound wave determines how much the light is refracted in the crystal. If the amplitude becomes very low, only the unrefracted beam is visible. The velocity with which the interference pattern moves depends on the derivative of the modulation function. For example, a triangle function leads to linear expansion/contraction, while the non-differentiable block function leads to sudden jumps in the positions of the fringes.

Now we observe the effect on the interference pattern if we change the phase difference. For a phase difference of  $0^\circ$  both patterns contract/expand at the same times. For  $180^\circ$ , they move exactly opposite to one another. For  $90^\circ$ , the maximum velocity of one pattern is reached when the other pattern is momentarily at rest.

A circle can be drawn at fringe (1, 1), if the phase difference between the AOMs is exactly  $90^\circ$  and they both have the same amplitude. To draw an ellipse, both AOMs can be given different amplitudes. If one wants to draw a circle with e.g. (1, 2), the amplitudes have to be adjusted accordingly by the same factor.

#### relation to polarization

The phase shifting box has a similar effect to that of a waveplate, since the phase of one wave is shifted relative to another.

If one frequency is an exact multiple of the other, Lissajous figures can be created, for example an eight figure for 2:1. If the ratio between the frequencies is not a natural number, then

the Lissajous figure is not closed, but open. If the ratio between the frequencies is not a natural number, then the Lissajous figure is not closed, but open

Now we again want to draw a circle.

The difference between a laser and light from e.g. a light bulb is that the light from the bulb is not linearly polarized, doesn't have a uniform frequency and is not coherent, i.e. doesn't have the same phase for all photons. An ideal laser would have these properties. coherence length/time?

## 4 Results

## 5 Analysis

## 6 Summary

## 7 Discussion



## References

- [1] lab course script for F85  
<https://www.physi.uni-heidelberg.de/Einrichtungen/FP/anleitungen/F85.pdf>